

THERMAL CYCLING CONDITION

STEPS	TEMP (°C)	TIME	QUANTITATION	CYCLE
Enzyme Activation	95	1 min	Off	1
PCR and Quantitation	95	5 sec	Off	40
	56	30 sec	On	

Table 4:- Thermal Cycling Condition

STANDARD CURVE CONTROL PARAMETER

Control Parameters	Valid Value
Slope	-3.00 / 3.74
PCR Efficiency	85% / 115%
R square (R²)	0.95-1

Table 5:- Curve Control Parameter

SAMPLE TYPE

- Cytomegalovirus (CMV) Quantitative Real-Time PCR Detection Kit developed by NeoDx Biotech Labs Pvt Ltd is a robust and convenient molecular diagnostic device for the Quantitation of CMV DNA in human Cerebrospinal Fluid, Plasma, Serum or Whole Blood DNA extracts.
- The kit can be used for Quantitation of CMV viral DNA if coupled with either NucleoDx Blood DNA or any other CE-IVD, EU or FDA approved DNA extraction kits.
- While using other DNA extraction kits, elute the DNA in 30-40 µL of elution buffer for accurate results.
- Always use freshly extracted patient DNA, or the DNA samples stored following proper guidelines, to get the optimal results.

RESULT VALIDATION

CMV Genes	CMV Assay Result
CMV DNA > 12 IU/mL	The result is within the acceptance range. The detection probability if CMV DNA is >95%. The Positive test result is statistically ensured.
CMV DNA < 12 IU/mL	The result is outside the acceptance range. The reproducibility of the positive result is not assured.
CMV DNA Negative	No CMV DNA was detected.

Table 6:- Result Validation

DATA ANALYSIS

- After completion of the run, analyze the data as per the Real-Time PCR manufacturer instructions.
- Analysis should be performed separately for each target using a manual threshold settings.
- In case internal control has not worked for a sample re-do the test with 2 or more dilutions.
- Negative results do not exclude infection and should not be used as the sole basis for the treatment.
- The quantitation standards are defined as IU/µL. The following equation has to be applied to convert the values determined using the standard curve into IU/ mL of sample material.

$$\text{Result (IU/mL)} = \frac{\text{Elution Volume (}\mu\text{L)} \times \text{Result (IU/}\mu\text{L)}}{\text{Sample Volume (mL)}}$$

Sample Volume is the volume of sample used for DNA extraction.

QUANTITATIVE STANDARD IN IU/µL

Standards	International Units (IU/ µL)
Quantitation Standard 1 (QS 1)	120,000
Quantitation Standard 2 (QS 2)	12,000
Quantitation Standard 3 (QS 3)	1200
Quantitation Standard 4 (QS 4)	120
Quantitation Standard 5 (QS 5)	12

Table 7:- Quantitative Standards

WARNING

Positive control showed no amplification

Inappropriate storage of reagents

- Store the reagents at recommended temperature for their optimal performance.
- Avoid repetitive freezing and thawing.
- Check the expiry of reagents.

Negative controls are positive

Causes - Cross-Contamination

- Follow good laboratory practices to avoid contamination issues.
- Use a new batch of reagents and repeat the experiment.

Abnormal plot and/or low ΔRn values in amplification curve

The baseline was set improperly (some samples have CT values lower than the baseline value)

- Switch from manual to automatic baseline, or move the baseline stop value to a lower CT (2 cycles before the amplification curve for the sample crosses the threshold).

An amplification signal is detected in the early cycles

Dilute the sample to increase the CT value.

REFERENCES

Optimization of Quantitative Detection of Cytomegalovirus DNA in Plasma by Real-Time PCR by Michael Boeckh , MeeiLi Huang , James Ferrenberg , Terry Stevens-Ayers , Laurence Stensland , W. Garrett Nichols , and Lawrence Corey
<https://doi.org/10.2174/187152611797636703>